

Neb-u-Tyke Bella Butterfly

This interactive, child-friendly design provides an ideal solution for effective aerosol treatment for children.

Designed with an attractive butterfly shape and colorful stickers to personalize, Bella Butterfly will keep child relaxed and entertained during aerosol treatment. The John Bunn Neb-u-Tyke Bella Butterfly is the latest addition to the Graham-Field line of quality nebulizer compressors that offer innovation, performance, and value for both patients and providers.


Features:

- Reliable, piston-driven compressor
- Child-friendly design complete with stickers to personalize during aerosol treatment
- Convenient and compact design for easy transport and storage
- Comes complete with Nebulizer Accessory Kit
- Five year limited warranty
- Designed to meet the requirements of HCPCS Code: E0570

The John Bunn Neb-u-Tyke Bella Butterfly Nebulizer comes complete with:

- 1 Neb-u-Tyke Bella Butterfly Nebulizer Compressor
- 1 User Manual
- 1 Nebulizer Accessory Kit complete with
- 1 Reusable nebulizer (single patient use only)
- 1 Pediatric mask
- 1 Mouthpiece
- 1 PVC tubing
- 5 Spare filters

JB0112-071	Neb-u-Tyke Bella Butterfly Nebulizer	1 ea
JB0112-071F	Replacement filters	5/pk



Specifications

Power	120VAC, 60Hz	Weight	2.88 lb (1.31 kg)
Operating Temperature Range	50°F TO 104°F	Max Pressure (Compressor)	42 psi
Operating Humidity Range	10 to 95% RH	Max Flow (Compressor)	Approx 8 LPM
Storage Temperature Range	4°F TO 140°F	Particle Size (Nebulizer)	≤5 microns
Storage Humidity Range	10 to 95% RH	Capacity (Nebulizer)5 ml
Dimensions.....	6.34" x 9.1" x 3.7" (161 mm x 231 mm x 93 mm)	Method Of Nebulization.....	Piston Driven
		Sound Level.....	under 60 dBA

Specifications are subject to change without notice.

FDA (USA) restricts the sale of this device to or by the order of a physician.