

Neb-u-Tyke Neb-a-Doodle™

This interactive, child-friendly design provides a solution for effective aerosol treatment for children. Designed with an integrated drawing board, the Neb-a-Doodle will keep child relaxed and entertained during aerosol treatment and can be used over and over again. The John Bunn Neb-u-Tyke Neb-a-Doodle is the latest addition to the Graham-Field line of quality nebulizer compressors that offer innovation, performance, and value for both patients and providers.

Features:

- Reliable, piston-driven compressor
- Interactive, drawing board for child to draw, doodle and write during aerosol treatment
- Convenient and compact design for easy transport and storage
- Comes complete with Nebulizer Accessory Kit
- Five year limited warranty
- Designed to meet the requirements of HCPCS Code: E0570

The John Bunn Neb-u-Tyke Neb-a-Doodle Nebulizer comes complete with:

- 1 Neb-u-Tyke Neb-a-Doodle Nebulizer Compressor
- 1 User Manual
- 1 Nebulizer Accessory Kit complete with
 - 1 Reusable nebulizer (single patient use only)
 - 1 Adult mask
 - 1 Pediatric mask
 - 1 Mouthpiece
 - 1 PVC tubing
 - 1 "T" adapter
 - 3 Spare filters
 - 1 Reservoir tube



- JB0112-070** Neb-u-Tyke Neb-a-Doodle
- JB77400** Neb-u-Tyke Neb-a-Doodle & Case
- JB76240** Case only



Specifications

Power	115VAC, 60Hz	Weight	4 lbs
Power Consumption	≤ 60 watts	Max Pressure (Compressor)	32.5 psi
Operating Temperature Range	41°F TO 104°F	Max Flow (Compressor)	Approx 12 LPM
Operating Humidity Range	30 to 95% RH	Particle Size (Nebulizer)	≤5 microns
Storage Temperature Range	40°F TO 158°F	Capacity (Nebulizer)	10 ml
Storage Humidity Range	10 to 95% RH	Method Of Nebulization	Piston Driven
Dimensions	10" × 6.5" × 4.2"	Sound Level	under 57 dBA
	(255 mm x 166m x 108 mm)		

*Specifications are subject to change without notice.
 FDA (USA) restricts the sale of this device to or by
 the order of a physician.*