



Feather® Sterile Scalpels

- Sterile, disposable scalpel features heat-treated stainless steel blade securely integrated with rigid, chemically resistant plastic handle.
- Individually wrapped in hermetically-sealed aluminum foil.
- Blades are protected with snap-off cover.
- Scalpels available with choice of blade size.
- All sizes packed 20/bx.

- 2975#10
- 2975#11
- 2975#12
- 2975#15
- 2975#20
- 2975#21
- 2975#22
- 2975#23
- 2975#24



Feather Sterile Scalpels

Stainless Steel Blades

#10	#11	#12
#15	#20	#21
#22	#23	#24