

- Leg and forearm sections can be adjusted independently for a custom fit
- Vinyl-coated arm cuffs can be shaped to the user's arm for extra comfort
- Contoured vinyl handgrip is durable and comfortable
- Telescoping parts have internal bushings and external lock nuts to silence any noise
- Extra-large, non-skid crutch tips provide added traction and stability
- Constructed with high strength aluminum tubing
- Easy push-button height adjustment with locking collet for added security
- Grip-to-floor height adjusts in 1" increments
- Arm cuff to grip height adjusts in 1" increments
- Includes laminated color operating instructions
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0110
- Three (3) year limited warranty

6341A	Large	1 pr
6342A	Medium	1 pr
6343A	Small	1 pr

