



- Made of lightweight, anodized aluminum
- Reinforced center tube and footpiece provide added strength
- Push button adjustment makes changing height easy
- Height adjustment in 1" increments
- Handgrip position is readily adjusted with wing nuts
- Underarm and handgrip pads are comfortable and durable
- Includes laminated color operating instructions
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0114
- Three (3) year limited warranty

- 3610LF-8** Adult 8 pr/cs
3611LF-8 Tall 8 pr/cs
3612LF-8 Youth 8 pr/cs
3613C-8 4 Adult, 2 Tall, 2 Youth 8 pr/cs

	Armpad to Grip Height Range	Armpad to Floor Height Range	Recommended Patient Height Range
3610LF-8	15" – 19"	45" – 53"	5' 2" – 5' 10"
3611LF-8	16" – 20"	53" – 61 1/2"	5' 10" – 6' 6"
3612LF-8	9" – 13"	37" – 46"	4' 6" – 5' 2"

