

Economy Overbed Table, Non-Tilt

- Easy to assemble
- Height adjusts infinitely between 28" and 44"
- Twin-wheel swivel casters provide smooth-rolling action. Two locking and two non-locking casters
- Walnut, wood-grain laminated ¾" particle-board core with plastic edge molding
- Chrome plated "H" base and column for aesthetics and durability
- Designed to meet the requirements of HCPCS code: E0274
- Supports up to 30 lbs. of dispersed weight.
- Warranty: One Year

GF8902 1 ea

Specifications

Table Height.....	28 - 44"
Table Width.....	30"
Table Depth.....	15"
Base Width.....	26.2"
Base Depth.....	15.25"
Weight.....	22 lbs
Dispersed Weight Capacity.....	30 lbs

