

**Fast Alert Advanced Patient Alarm**

- Alarm notifies staff or caregiver when user gets out of bed or chair
- Pressure sensitive pad with 1 year warranty connects to alarm
- Alarm easily secures with self contained clip onto areas such as wheelchairs
- Convenient dual volume settings: 70dB and 90dB
- Delay function allows for 0 or 2 second delay to reduce likelihood of false alarms with restless users
- Reset button on side of monitor silences and resets alarm
- Remote alarm option allows the caregiver to silence the local alarm and sound the alarm remotely through the nurse call system (Adapter Cable and Splitter sold separately)
- Requires one 9V battery (not included)
- Comes with protective boot cover to prevent damage if dropped
- Off switch and Low battery indicator
- 1 Year Limited Warranty
- Pad size:  
 Chair - 10" x 15"  
 Bed - 10" x 29"

<b>GF13702</b>	Advanced Patient Alarm only	1 ea
<b>GF13702B</b>	Advanced Patient Alarm with Bed Pad	1 set
<b>GF13702C</b>	Advanced Patient Alarm with Chair Pad	1 set

*Replacement Parts*

<b>106375</b>	Pressure Sensor Bed Pad (1 Year Warranty)	1 ea
<b>106380</b>	Pressure Sensor Chair Pad (1 Year Warranty)	1 ea
<b>GF1370-BCLIP</b>	Bed Clip	1 ea
<b>GF1370-COVER</b>	Replacement Protective Boot Cover	1 ea
<b>106256</b>	Nurse Call Splitter	1 ea
<b>105681</b>	Adapter cable for Alarm	1 ea
<b>105682</b>	AC Adapter for Alarm	1 ea

