

**Floor Stand**

- Provides freestanding capability for Versa-Helper Trapeze 2800A and 2800GA.
- Assisted weight capacity of 250 lbs.
- Individually-leveled floor bars provide firm, secure base.
- Folds flat for easy storage and handling.
- When used with Versa-Helper Trapeze, designed to meet the requirements of HCPCS code: E0940.

<b>2840A</b>	Chrome-plated	1 ea
<b>2840GA</b>	Gray	1 ea

**Specification**

Base Length .....	40"
Base Width .....	29.2"
Width Between Legs.....	27"
Height .....	39"

