




Clinical Care 600S Stretcher Pad

The Simmons Clinical Care stretcher pad is the best choice when selecting a non-powered support surface. Utilizing a tri-laminate design in combination with high resiliency foam and contour cut surface modification, the Simmons Clinical Care 600S immerses the user while providing durable support.

- 5" height for maximum support
- Contour cut grid allows for better patient immersion and greater airflow
- Better pressure redistribution properties than standard pad offerings
- Tri-laminate design for support and durability
- Four-way stretch cover material is fluid resistant and anti-microbial
- Available in many sizes and configurations to match your stretcher model
- 1 Year Warranty



Not made
with natural
rubber latex



SIMMONS Clinical Care products manufactured by GF Health Products, Inc. under license from the trademark owner, Dreamwell, Ltd., a subsidiary of Simmons Bedding Company, Atlanta, Georgia.

SIMMONS is a trademark of Dreamwell, Ltd. a subsidiary of Simmons Bedding Company.

© 2014 Dreamwell, Ltd. All rights reserved.



www.grahamfield.com

© 2014 GF Health Products, Inc.
Atlanta, GA 30360 All rights reserved.
GF1300126RevB14