



Lightweight Single Head Stethoscope

The Labtron Lightweight Single Head Stethoscope is ideal for general assessment and blood pressure measurement. The slim, lightweight anodized aluminum chestpiece easily slides beneath a blood pressure cuff. Labtron stethoscope headsets are anatomically adjustable to ensure a comfortable acoustically sealed aural fit.



Features:

- Complete unit weighs only 2.5 oz (70g)
- 1 3/4" (4.5cm) diameter diaphragm with non-chill diaphragm ring
- Adjustable anodized aluminum binaurals
- Flexible PVC tubing
- 22" (56.5cm) tubing; 31.5" (80cm) total length
- Full color retail packaging
- Includes stethoscope with diaphragm only chestpiece, adjustable binaurals, and (1) set of white firm threaded PVC eartips
- One (1) year limited warranty



| Item | Color | UOM |
|-----------|------------|------|
| 300DLX | Black | 1 ea |
| 300DLX-BL | Blue | 1 ea |
| 300DLX-GL | Gold | 1 ea |
| 300DLX-GY | Grey | 1 ea |
| 300DLX-R | Red | 1 ea |
| 300DLX-Y | Yellow | 1 ea |
| 300DLX-LB | Light Blue | 1 ea |
| 300DLX-TL | Teal | 1 ea |



**This product is not made with natural rubber latex*



www.grahamfield.com

© 2016 GF Health Products, Inc.
Atlanta, GA 30360 All rights reserved.
GF1600019RevA16