



A Graham-Field Family Brand

S300 SERIES

- Low acuity mattress is first line of defense against patients “at-risk” for pressure ulcers
- “Flat-Foam” top surface provides comfort and pressure reduction
- Tri-Laminate construction with soft heel-zone
- Firm bolster with optional 3" raised bolster
- Meets 16 CFR 1633 Flammability Standard
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED

Warranty:

Mattress: Three (3) year limited warranty

Cover: One (1) year limited warranty



LS
*Latex Statement

**This product is not made with natural rubber latex*

Dimensions	Flush Bolster	Raised Bolster
76" X 35"	S300B-3576	S300RB-3576
80" X 35"	S300B-3580	S300RB-3580
84" X 35"	S300B-3584	S300RB-3584
76" X 42"	S300B-4276	S300RB-4276
80" X 42"	S300B-4280	S300RB-4280
84" X 42"	S300B-4284	S300RB-4284



SIMMONS Clinical Care products manufactured by GF Health Products, Inc. under license from the trademark owner, Dreamwell, Ltd., a subsidiary of Simmons Bedding Company, Atlanta, Georgia. SIMMONS is a trademark of Dreamwell, Ltd. a subsidiary of Simmons Bedding Company.

© 2016 Dreamwell, Ltd. All rights reserved.



www.grahamfield.com



© 2016 GF Health Products, Inc.
Atlanta, GA 30360 All rights reserved.
GF1300001RevB16