





# S500 SERIES

A Graham-Field Family Brand

- Preferred combination of pressure redistribution and support
- Modified convoluted topper enhances immersion characteristics and increases airflow
- Tri-laminate construction with soft heel-zone
- Firm bolster with optional 3" raised bolster
- Fluid-resistant, low friction, anti-microbial top cover with vinyl bottom cover and full zipper
- Designed to meet the requirements of HCPCS Code E0184
- 400 lb (450 lb for 42" wide models) maximum weight capacity, EVENLY DISTRIBUTED
- Meets 16 CFR 1633 Flammability Standard

**Warranty:**

Mattress: Five (5) year limited warranty  
 Cover: One (1) year limited warranty



*\*This product is not made with natural rubber latex*

Dimensions	Flush Bolster	Raised Bolster
76" X 35"	S500B-3576	S500RB-3576
80" X 35"	S500B-3580	S500RB-3580
84" X 35"	S500B-3584	S500RB-3584
76" X 42"	S500B-4276	S500RB-4276
80" X 42"	S500B-4280	S500RB-4280
84" X 42"	S500B-4284	S500RB-4284



*SIMMONS Clinical Care products manufactured by GF Health Products, Inc. under license from the trademark owner, Dreamwell, Ltd., a subsidiary of Simmons Bedding Company, Atlanta, Georgia.  
 SIMMONS is a trademark of Dreamwell, Ltd. a subsidiary of Simmons Bedding Company.  
 © 2016 Dreamwell, Ltd. All rights reserved.*



[www.grahamfield.com](http://www.grahamfield.com)



© 2016 GF Health Products, Inc.  
 Atlanta, GA 30360 All rights reserved.  
 GF1200156RevB16