


### Tuning Fork – Adjustable Frequency

- All aluminum tuning fork, with adjustable weight scale in seven frequencies: C128, D144, E160, F170.7, G192, A213.3 and B240 (weights removed).
- Scale letters are stamped on stems.

1325    Tines 3/8"    Length 8"    1 ea



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**[www.grahamfield.com](http://www.grahamfield.com)**