

### Standard Finger Ring Cutter

- Specially designed lever slips under ring to permit rotating thumb screw to cut ring.
- Knurled straight handle with 7/8" rotary jeweler's saw.

**1372-3**      5 3/4"      1 ea

**1372-4**      Ring Rotary Blade Only      1 ea



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**[www.grahamfield.com](http://www.grahamfield.com)**