

## Features:

- Offset handle places the user's weight over the base for maximum balance and control
- Low profile, 4-point base further enhances stability
- Rounded edge base plate helps prevent scraping and bruising of ankles
- Attractive finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Height adjusts from 29" - 38" in 1" increments
- Tool-free left to right hand adjustment is made easily
- Choice of silver or bronze finish
- Includes laminated color operating instructions
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0105
- Three (3) year limited warranty

### Large Base (12" x 8")

<b>6140A</b>	Aluminum	4/cs
<b>6140BZ</b>	Bronze	4/cs
<b>6140BK</b>	Black	4/cs

### Small Base (8" x 6")

<b>6141A</b>	Aluminum	4/cs
<b>6141BZ</b>	Bronze	4/cs
<b>6141BK</b>	Black	4/cs

