

Silver Collection Low Profile Quad Canes with Ortho-Ease® Grip

Ortho-Ease® Grip

- Grip angle is designed to adjust to comfortably fit the user's hand and wrist angle and to minimize stress on the joints
- Ortho-Ease® grip is ideal for arthritics and for those needing a wider grip for greater support
- Comes complete with convenient wrist strap

Features:

- Offset handle places the user's weight over the base for maximum balance and control
- Low profile, 4-point base further enhances stability
- Rounded edge base plate helps prevent scraping and bruising of ankles
- Attractive finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Height adjusts from 29" 38" in 1" increments
- Tool-free left to right hand adjustment is made easily
- Includes laminated color operating instructions

1 ea

- Designed to meet the requirements of HCPCS code: E0105
- Maximum weight capacity 300 lb EVENLY DISTRIBUTED
- Three (3) year limited warranty

Large Base (12" x 8")
6120A Aluminum 4/cs
Small Base (8" x 6")
6121A Aluminum 4/cs

6121A-1 Aluminum



