

LUMEX®

A Graham-Field® Family Brand

Adjustable Offset Canes - Nitrile Grip



Features:

- Offset handle places the user's weight directly over the base for maximum balance and control
- Attractive finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Available in a variety of finishes
- Includes laminated color operating instructions
- Soft, warm and contoured for maximum comfort
- Sculpted grip reduces the chance of slippage, even when wet
- Comes complete with convenient wrist strap
- Designed to meet the requirements of HCPCS code: E0100
- Three (3) year limited warranty
- 250 lb maximum weight capacity - EVENLY DISTRIBUTED



5941A	Bronze	Standard Length	(31"-39")	6/cs
6333A	Bronze	Extra Long	(34.5"-42.5")	1 ea
5942A	Black	Standard Length	(31"-39")	6/cs
5943A	Silver	Standard Length	(31"-39")	6/cs
5944A	Multi-Pack	Standard Length	(31"-39")	6/cs
	(2) Bronze, (2) Black, (2) Silver			

graham
field™

www.grahamfield.com

© 2016 GF Health Products, Inc.
Atlanta, GA 30360 All rights reserved.
GF1500070RevB16