

The Advantage Recliner comes with a higher back and removable full support headrest. The headrest includes an additional neck support that adjusts easily to fit every need. The durable hammertone frame is designed to support up to 300 lb and in the bariatric version up to 450 lb users.

Reclined position



Standard Features

- Available in 16", 18" and 20" seat widths, with 17" seat depth and in 22" seat width with 19" seat depth
- Features manual full reclining back range from 90-170 degrees
- Removable full support headrest
- Detachable padded desk or full arms
- Hammertone finish provides durability and style
- Embossed black leatherette type vinyl upholstery for easy cleaning
- Anti-tippers standard
- Comes with spreader metallic back bar which works as anti-fold device
- Industry compatible tool-free adjustable elevating front rigging
- Standard chart pocket on back
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED for 16, 18 and 20" models
- Maximum weight capacity 450 lb - EVENLY DISTRIBUTED for 22" models
- 16", 18", 20" models meet Medicare code K0001/ E1226 (Recliner Back) + E0966 (Head Rest Extension) + E0971 (Anti-tipping device) + K0195 (Elevating Legrests)
- 22" models meet Medicare code K0007/E1226 (Recliner Back + E0966 (Head Rest Extension) + E0971 (Anti-tipping device) + K0195 (Elevating Legrests) + E2201

Warranty

Five (5) year limited warranty





Headrest



Chart pocket and spreader bar



Anti tipper

Advantage Recliner		Dimensions and Weights									
		Seat			Back	Overall			Wt.	Wt. Cap	
3K010130	Desk Arm, Elevating Legrest	18"	17"	19 1/2"	23"	34"	51"	26 1/2"	11 1/2"	53 lb	300 lb
3K010150	Full Arm, Elevating Legrest	18"	17"	19 1/2"	23"	34"	51"	26 1/2"	11 1/2"	53 lb	300 lb
3K010230	Desk Arm, Elevating Legrest	16"	17"	19 1/2"	23"	34"	51"	28 1/2"	11 1/2"	53 lb	300 lb
3K010250	Full Arm, Elevating Legrest	16"	17"	19 1/2"	23"	34"	51"	28 1/2"	11 1/2"	53 lb	300 lb
3K010330	Desk Arm, Elevating Legrest	20"	17"	19 1/2"	23"	34"	51"	30 1/2"	11 1/2"	53 lb	300 lb
3K010350	Full Arm, Elevating Legrest	20"	17"	19 1/2"	23"	34"	51"	30 1/2"	11 1/2"	53 lb	300 lb
3K020430	Desk Arm, Elevating Legrest	22"	19"	19 1/2"	23"	34"	51"	32 1/2"	11 1/2"	57 lb	450 lb
3K020450	Full Arm, Elevating Legrest	22"	19"	19 1/2"	23"	34"	51"	32 1/2"	11 1/2"	57 lb	450 lb

Advantage Recliner – Common Accessories & Parts		HCPCS Code
90763062	Wheel lock extension, pair	E0961
90763064	Heel loop, pair	E0951
92133040	Positioning belt 37" (Airline style) (Metal buckle)	E0978
90133125	Positioning belt 48" long (Airline style) (Metal buckle)	E0978
300-505	Positioning belt (Auto-style) push button	E0978
90763730	Anti-tippers, pair for 16x17/18x17/20x17 (300 lb)	E0971
90763067	Anti-tippers, pair for 22x19 (450 lb)	E0971
92523850	Adjustable chromed anti-fold bar (lower frame) chromed	
907637E1	Hammertone back spreader bar 16" W (reclining back)	
907637E2	Hammertone back spreader bar 18" W (reclining back)	
907637E3	Hammertone back spreader bar 20" W (reclining back)	
907637E4	Hammertone back spreader bar 22" W (reclining back)	
Complete Headrest Extension is composed by: Posts, cushion and headrest upholstery (ordered separately)		
90763720	Headrest posts, pair	
90763735	Headrest cushion, black vinyl	
907637D1	Headrest upholstery black 16" W (includes hardware)	
907637D2	Headrest upholstery black 18" W (includes hardware)	
907637D3	Headrest upholstery black 20" W (includes hardware)	
907637D4	Headrest upholstery black 22" W (includes hardware)	