

Grafco® Unlabeled Flint Glass Applicator and Dressing Jars are excellent for storing and dispensing a variety of supplies, making them ideal for hospitals and clinics.



Features:

- Transparent flint glass construction provides easy visual of contents inside
- Each jar includes removable overlapping polished aluminum cover
- Reusable
- Standard jars not autoclavable
- Pyrex® jars are sterilizable and autoclavable
- One (1) Year Limited Warranty

Item	Description	Dimensions	Case Weight	UOM
3459P	Pyrex Applicator Jar	6 5/8" x 3"	.7 lb per jar	1 ea
3460P	Pyrex Dressing Jar	3" x 3"	3.9 lb per case	12/cs
3461P	Pyrex Dressing Jar	4" x 4"	4.4 lb per case	6/cs
3462P	Pyrex Dressing Jar	5" x 5"	5.1 lb per case	6/cs
3463P	Pyrex Dressing Jar	6" x 6"	18.4 lb per case	6/cs
3459	Standard Applicator Jar	6 5/8" x 3"	.7 lb per jar	1 ea
3460	Standard Dressing Jar	3" x 3"	3.9 lb per case	12/cs
3461	Standard Dressing Jar	4" x 4"	4.4 lb per case	6/cs
3462	Standard Dressing Jar	5" x 5"	5.1 lb per case	6/cs
3463	Standard Dressing Jar	6" x 6"	18.4 lb per case	6/cs



**This product is not made with natural rubber latex*