


The Lumex Bedside Assist Rail offers assistance to the patient/resident for getting in and out of a home style bed.



- Can be used on either side of the bed for assistance with getting in and out of home style bed
- Strong, steel construction with attractive black finish makes it durable and look good in the home
- Cushioned handle that is adjustable in height assists with getting in and out of bed
- Easily slides in between mattress and bed frame with added feature of anti-slip foam on base bar and securing strap to secure to bed frame
- Convenient, black pouch included for storage of personal items
- Tool-free assembly
- One (1) year limited warranty



GF6900 Black 1 ea

Specifications

| | |
|-----------------------------------------------------------------------|---------------|
| Height adjustment from mattress base to top of handle..... | 17.5" – 21.5" |
| Length of rail base tubing support under mattress | 36" |
| Overall width of handle..... | 19.5" |
| Product weight | 6.7 lb |
| Weight capacity:..... | 300 lb |
| (product is not designed to support the full weight of an individual) | |