

The Walkabout ConTour Imperial Hemi is the hemi solution for the bariatric market! It has the same great features as the ConTour Imperial and folds quickly and easily for storage and travel. Features include a padded seat, ergonomic hand grips, easy-to-operate hand brakes, and a contoured, padded backrest.



## Features:

- 20.5" seat to floor height
- Epoxy coated aluminum frame
- Locking loop brakes
- Removable pouch under seat for storage of personal items
- Adjustable handle height
- Heavy-duty 8" wheels for indoor and outdoor use on even terrain
- 450 lb maximum weight capacity - EVENLY DISTRIBUTED
- Weighs 17.5 lb
- Designed to meet the requirements of HCPCS codes: E0149 + E0156
- Includes laminated color operating instructions
- Five (5) year limited warranty



**RJ4402R** Burgundy 1 ea

**Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.**

## Specifications

Width Inside Hand Grips .....	21.5"
Overall Width .....	28.75"
Overall Depth.....	26"
Seat Width.....	19.5"
Seat Depth .....	13.8"
Seat Height.....	20.5"
Caster Diameter .....	8"
Minimum Handle Height.....	31"
Maximum Handle Height.....	34.5"
Weight Capacity .....	450 lb
	EVENLY DISTRIBUTED