

**FR565TG PREFERRED CARE RECLINER
TILT-IN-SPACE: REPLACEMENT PARTS LIST**

TO ORDER: Tel: 1.770.368.4700 • Fax: 1.770.368.2386 • www.grahamfield.com

Returns: Parts may not be returned without prior authorization from GF Health Products, Inc. Custom parts cannot be returned for credit. Please refer to GF Health Products, Inc.

Terms and Conditions for further information.



ITEM	ITEM NUMBER	DESCRIPTION	UNIT OF MEASURE
1	5654A825G	DROP ARM ASSEMBLY LEFT, FOR FR565DG / FR565TG	1 EA
2	5654A815G	DROP ARM ASSEMBLY RIGHT, FOR FR565DG / FR565TG	1 EA
3	56575U(COLOR CODE)	ARMREST, UPHOL FOR FR565DG AND FR565TG	1 EA
4	56575U(COLOR CODE)FR	ARMREST, UPHOL FOR FR565DG AND FR565TG	1 EA
5	56548U(COLOR CODE)FR	SIDE PANEL, LEFT UPHOL, FR FOR FR565DG / FR565TG	1 EA
6	56549U(COLOR CODE)FR	SIDE PANEL, RIGHT UPHOL, FOR FR565DG / FR565TG	1 EA
7	56558J(COLOR CODE)FR	ARMREST AND SIDE PANEL, LEFT, UPHOLSTERED FOR FR565WG SERIES	1 EA
8	56559J(COLOR CODE)FR	ARMREST AND SIDE PANEL, RIGHT, UPHOLSTERED FOR FR565WG SERIES	1 EA
9	56558U(COLOR CODE)FR	ARMREST AND SIDE PANEL, LEFT, UPHOLSTERED FOR FR565G SERIES	1 EA
10	56559U(COLOR CODE)FR	ARMREST AND SIDE PANEL, RIGHT, UPHOLSTERED FOR FR565G SERIES	1 EA
11	56512U(COLOR CODE)	BACK, UPHOLSTERED FOR 565G	1 EA
12	56512U(COLOR CODE)FR	BACK, UPHOLSTERED FR FOR FR565G	1 EA
13	5651A300G	SEAT AND LEGREST, UPHOLSTERED FOR 565G	1 EA
14	56511U(COLOR CODE)FR	SEAT AND LEGREST, UPHOLSTERED FR FOR 565G	1 EA



ITEM	ITEM NUMBER	DESCRIPTION	UNIT OF MEASURE
15	WC150006	CASTER, FRONT FOR FR565 SERIES, NOT FR565WG	1 EA
16	WW170007	CASTER, REAR LOCKING FOR FR565DG / FR565TG SERIES, NOT FR565WG	1 EA
17	WB190002	FOOT PEDAL, FOR FR565TG	1 EA
18	WB190004	WEBBING, BACK FOR FR565 SERIES	1 EA
19	5655A001	CABLE ASSEMBLY, BACK FOR FR565TG (not shown)	1 EA
20	HP327220	HARDWARE, CLEVIS PIN FOR GAS SPRING AND FR565 SERIES BACK	1 EA
21	5792A003	HARDWARE, CLEVIS PIN W COTTER CLIP 577 GAS SPRING AND 565 SERIES BACK	10 EA / ST
22	5770A073	HARDWARE, GAS SPRING BRACKET	1 EA
23	EI000000	HARDWARE, GAS SPRING FOR FR565TG	1 EA
24	EI000001	HARDWARE, GAS SPRING FOR FR577RG / FR587W / FR565TG	1 EA
25	VM000003	HARDWARE, TORQUE LIMITER FOR FR565TG (not shown)	1 EA
26	LRMP0001	PAINT TOUCH-UP LUMEX GRAY (not shown)	1 EA



1.770.368.4700

Information contained herein is subject to change.
 The most current and complete product information can be found on our website.
www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
 Graham-Field and Lumex are trademarks of GF Health Products, Inc.
 GF Health Products, Inc. is an ISO 13485 Certified Company.