


## PVC Reclining Shower Chair

- Constructed of the highest quality health care grade PVC Dual reinforcement into the caster fittings adds additional support
- Quick drying mesh available in multiple colors.
- Threaded stem casters for additional safety.
- Ergonomic curved rails provide support without sharp edges and reclining chairs have anti-slip handgrips for added safety.
- One year limited warranty on PVC frame; six month limited warranty on all other materials.
- Deluxe elongated extra heavy-duty, contoured open front soft seat.
- Padded lumbar support and head section.
- Dual drop arms for patient transfer can be secured with red safety straps

89330 1 ea  
89440 1 ea  
89550 1 ea



89330



89440

Item#	Internal Width	External Width	Height	Depth	Seat Height	Back Width Exterior	Exterior Length Upright (leg down)	Exterior Length Reclined (leg up)	Reclined Back	Weight Capacity (lbs.)	Caster
89330	20"	24"	48"	18"	21"	16"	43"	58"	29"	325	heavy duty
89440	26"	30"	50"	22"	22"	23¼"	41"	62"	27"	600	heavy duty
89550	30"	34"	50"	22"	22"	27¼"	41"	62"	27"	700	heavy duty



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

20170 Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**www.grahamfield.com**

GF0700129RevB07