

The Comfort Care 600C stretcher pad is the best choice when selecting a non-powered support surface. Utilizing a tri-laminate design in combination with high resiliency foam and contour cut surface modification, the Comfort Care 600C immerses the user while providing durable support.

- 5" height for maximum support
- Contour cut grid allows for better patient immersion and greater airflow
- Better pressure redistribution properties than standard pad offerings
- Tri-laminate design for support and durability
- Four-way stretch cover material is fluid resistant and anti-microbial
- Available in many sizes and configurations to match your stretcher model
- 1 Year Warranty



**This product is not made with natural rubber latex*