

Tri-Step Cane

The Lumex Tri-Step Cane is the ultimate cane solution designed to keep you on the go. The base provides three point contact with the ground which allows the Tri-Step to stand upright on its own and provides added stability. The contoured grip is comfortable for the hands. With its folding design and strap to secure, the cane is very convenient for storage. And to make it more eye appealing, the Tri-Step is available in 4 attractive colors.



Features:

- Three point contact base allows the cane to stand upright on its own
- Sides of base can be used to step on to lift cane if cane has fallen on the ground – no more bending over to pick up a fallen cane!
- Comfortable, contoured grip
- Lightweight aluminum, folding design makes it easy for storage
- Easily adjusts in height from 33" 37"
- Designed to support up to 300 lb
- Product weight: 1.1 lb
- Cane tip can be purchased separately to fit most standard canes with ¾" (19mm) tubing (Item # 6000T)
- Available in Black, Blue, Green and Pink
- Designed to meet the requirements of HCPCS E0100
- Three (3) year limited warranty

6000BK	Black	4 ea/cs
6000B	Blue	4 ea/cs
6000G	Green	4 ea/cs
6000P	Pink	4 ea/cs
6000T	Replacement Tip	1 ea





© 2017 GF Health Products, Inc. Atlanta, GA 30360 All rights reserved. GF1700049RevA17