



**This product is not made with natural rubber latex*

- Protection for higher risk patients with pressure ulcers and high-risk heels
- Provides pressure reduction and redistribution
- Firm side perimeter construction with optional 3" raised perimeter bolsters (RB) featuring ingress-egress cutout
- Three-Layer construction with heel slope
- Fluid-resistant, low friction, anti-microbial top cover with vinyl bottom cover and full zipper
- Designed to meet the requirements of HCPCS Code E0184
- Maximum Weight Capacity 450 lb (500 lb for 42" wide models)
- EVENLY DISTRIBUTED
- Five (5) year limited warranty on mattress,
One (1) year limited warranty on cover

Dimensions	With Firm Side Bolsters	With 3" Raised Firm Side Bolsters
35" x 76" x 6"	C600B-3576	C600RB-3576
35" x 80" x 6"	C600B-3580	C600RB-3580
35" x 84" x 6"	C600B-3584	C600RB-3584
42" x 76" x 7"	C600B-4276	C600RB-4276
42" x 80" x 7"	C600B-4280	C600RB-4280
42" x 84" x 7"	C600B-4284	C600RB-4284