

Designed to aid the transfer from wheelchair, bed, commode, or chair



Features:

- Solid wood construction
- Lacquer-coated finish for smooth transfer
- Tapered ends for easier transfer
- Available with one or two handles



Item	Size:	UOM:	Weight Capacity:
(2 Handles):			
5242-24A-1	24"L x 8"W	1 ea	440 lb
5242-30A-1	30"L x 8"W	1 ea	440 lb
5242-32A-1	32"L x 10"W	1 ea	735 lb
(1 Handle):			
5241-24A-1	24"L x 8"W	1 ea	440 lb
5241-30A-1	30"L x 8"W	1 ea	440 lb