




## Feather® Sterile MicroScalpels

- Lightweight stainless steel handle with non-slip finger grip combines with angles, ultra sharp blade.
- Individually wrapped to assure sterility; color-coded for identification.

2979#15	15° angle	5/bx
2979#30	30° angle	5/bx
2979#45	45° angle	5/bx



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**[www.grahamfield.com](http://www.grahamfield.com)**