

The Lumex Walkabout Four-Wheel Imperial (Bariatric) is ideal for larger individuals. Similar to the Walkabout Lite, this rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 20 lbs and supports a maximum weight capacity of 450 lb - **EVENLY DISTRIBUTED**. The Walkabout folds quickly and easily into a compact unit for storage and transport.



## Features

- Ergonomic hand grips
- Easy to operate locking and “quick stop” loop brakes (pull up to stop; push down to lock)
- Padded 20" wide seat for increased seating surface and added comfort
- Removable padded backbar
- Adjustable handle height to accommodate broad range of user heights
- Durable, heavy duty 8" wheel for indoor and outdoor use on even terrain
- Basket included that can be placed underneath seat for carrying of personal items
- Designed to meet the requirements of HCPCS codes: E0149 & E1056
- Includes easy to read, retail-friendly user manual
- Five (5) year limited warranty

And now with added rear reflectors that increase visibility at night!




### Straight Backbar

- |                |          |      |
|----------------|----------|------|
| <b>RJ4400K</b> | Black    | 1 ea |
| <b>RJ4400R</b> | Burgundy | 1 ea |

### Contoured Backbar

- |                |            |      |
|----------------|------------|------|
| <b>RJ4405B</b> | Royal Blue | 1 ea |
| <b>RJ4405K</b> | Black      | 1 ea |
| <b>RJ4405R</b> | Burgundy   | 1 ea |



Color options

### Specifications

Width Inside Handgrips .....	21.5"
Overall Width .....	28.7"
Overall Depth .....	26"
	<i>measured wheel to wheel</i>
Seat Width .....	19.5"
Seat Depth .....	13.8"
Seat Height .....	21.8"
Caster Diameter .....	8"
Minimum Handle Height .....	34"
Maximum Handle Height .....	39"
Weight Capacity .....	450 lb
	<b>EVENLY DISTRIBUTED</b>
Product Weight .....	20 lb

**Rollators are not intended to be propelled while seated.  
Do NOT use this rollator as a wheelchair or transport device.**