

Features:

- Sturdy 1" aluminum tubing provides maximum strength while remaining lightweight
- Dual-release folding mechanism enables user to fold walker sides independently
- Accommodates all leg attachments
- Designed to meet the requirements of HCPCS code: (without wheels) E0135 (with wheels) E0143
- One (1) year limited warranty



Without wheels

716070A-2	Adult	2 ea/cs
716070A-4	Adult	4 ea/cs
716170A-4	Junior	4 ea/cs

With 5" front fixed wheels

716270A-2	Adult	Aluminum	2/cs
716270A-4	Adult	Aluminum	4/cs
716370A-2	Junior	Aluminum	2/cs
716370A-4	Junior	Aluminum	4/cs

Specifications:

	Adult	Junior
Width Inside Handgrips	18	16 ½"
Width Inside Rear Legs	21	19"
Overall Width	22	19 ½"
Depth Front to Rear	17	15"
Depth Folded	4	4"
Minimum Height	32	25 ½"
Maximum Height	39	34 ¼"
Weight Capacity - EVENLY DISTRIBUTED	300 lb	300 lb