

## Features:

- Flame-retardant\* white vinyl panels
- Easy to set up and to dismantle for storage
- Vinyl panels are washable and replaceable on snap out rods
- 1" (25mm) anodized aluminum frame
- 2" (51mm) twin wheel casters allow for easy movement
- One (1) year limited warranty

\*NFPA 701 Test Method 1-2010 Edition



Item	Description	UOM
4296W	Three section	1 ea
4296P	Replacement Panel Set (3)	1 set
4299C	Replacement Casters (4)	1 set

## Specifications

Each Panel: ..... 49" x 31"  
(1244mm x 787mm)  
Overall: ..... 70.25" x 79"  
(1784 mm x 2006mm)  
Weight: ..... 11.8 lbs (5.35 kg)