

Walkabout Steel Knock Down Rollator

The Lumex Walkabout Steel Knock Down Rollator offers a comfortable, traditionally designed, value priced solution to mobility for someone looking for an easy to use walking aid. The frame is made of a highly durable steel with 6" easy to maneuver casters. The Walkabout quickly and easily folds into a compact unit for storage and transport.



Features:

- Ergonomic hand grips make them easy to grip
- Adjustable handle height to accommodate broad range of user heights
- Durable construction steel frame
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Removable, hinged, padded backrest can be folded up and down
- 6" casters for ease of maneuverability
- Pouch included for storage of personal items
- Designed to meet the requirements of HCPCS codes: E0143 & E01056
- Three (3) year limited warranty



**Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.**

Item	Color	UOM
RJ5500R	Burgundy	1 ea
RJ5500K	Black	1 ea

Specifications:

Width Inside Handgrips	15.5" (394mm)
Overall Width	23.5" (597mm)
Overall Depth	22.5" (572mm)
Seat Width	14" (356mm)
Seat Depth	12" (305mm)
Seat Height	20" (508mm)
Caster Diameter	6" (152mm)
Minimum Handle Height	30" (762mm)
Maximum Handle Height	34" (864mm)
Product Weight	18 lb (8.2 kg)
Maximum Weight Capacity	325 lb

EVENLY DISTRIBUTED