



EJ777-3 Bariatric Transport Chair

TO ORDER:
Tel: 1.770.368.4700
Fax: 1.770.368.2386
www.grahamfield.com

Item	Part Number	Description	U/M
1	EJ76-BELT	Positioning Belt with Plastic 3-Prong Buckle	1 ea
2	EJ777-3ARM	Flip Desk Arm with Arm Pad, Side Panel, and Hardware	1 pr
3	EJ777-3BACK	Back Upholstery, 24" Wide Black Nylon with Hardware	1 ea
4	EJ777-3BK	Attendant Wheel Lock with Cable	1 pr
5	EJ777-3BRP	Upper Back Release Paddle	1 pr
6	EJ777-3CST	8" Front Caster with Fork and Hardware	1 pr
7	EJ777-3FAB	Front Arm Guide and Bracket	1 pr
8	EJ777-3FR	Front Rigging (Footrest)	1 pr
9	EJ777-3FTPLT	Black Composite Footplate	1 pr
10	EJ777-3RAB	Back Arm Guide and Bracket	1 pr
11	EJ777-3RWH	12" Rear Wheel	1 ea
12	EJ777-3S	Back Arm Sleeve	1 pr
13	EJ777-3SEAT	Seat Upholstery, 24" Wide Black Nylon with Hardware	1 ea
14	EJ777-3ANT	Anti-tippers (not shown) for only those EJ777-3 models with serial numbers after IOS291	1 pr



1.770.368.4700

Information contained herein is subject to change.
The most current and complete product information can be found on our website.
www.grahamfield.com

Follow us on

© 2015, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a trademark of GF Health Products, Inc.,
One Graham-Field Way, Atlanta GA 30340-3140

GF Health Products, Inc. is an ISO 13485 Certified Company.