



Features:

- Designed for users with limited strength or dexterity who are unable to grip the walker
- Adjustable strap comfortably secures arm in the forearm rest
- Padded forearm rest can be pivoted and angled to fit user comfortably and to assist in lifting walker
- Padded handgrip provides comfort
- Universal clamp fits most walkers
- Tool-free installation and adjustment allows quick fitting for most patient sizes and types
- Designed to meet the requirements of HCPCS code: E0154
- Three (3) year limited warranty



| Item | UOM |
|-------|------|
| 6132B | 1 pr |

*** It is recommended that the platform walker attachments be installed on both sides for better stability.**