

The Lumex Allura/Allura LX Rollators are designed to offer comfort and a sleek solution to mobility. The frame is made of a lightweight aluminum material that is both durable and makes mobility easier for any occasion. Designed with a flexible removable back, height adjustable ergonomic handles and being offered in a Tall version allow for many customization options for a broad range of users. The Allura LX is designed with a 26" seat height and 41.7" maximum handle height making it ideal for users taller than 6'.



Features:

- Attractive, sleek, lightweight aluminum frame design
- Flexible, padded back support for comfort
- Ergonomic handgrips
- Cross brace design allows for side-to-side folding for storage and transport and offers added stability to rollator frame
- Handle height easily adjusts with simple and secure push buttons
- Wide 18" padded seat with fold-up handle offers added comfort
- Taller version with 26" seat height and 41.7" maximum handle height is ideal for users taller than 6'0"
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Large 8" casters for easy maneuverability and indoor or outdoor use
- Removable carry pouch for storage of personal items
- Comes with cane holder on side

Item	Description	Color	UOM
LX5000A	Standard	Apple Red	1 ea
LX5000BK	Standard	Black	1 ea
LX5000LXA	Tall	Apple Red	1 ea
LX5000LXBK	Tall	Black	1 ea

	LX5000	LX5000LX
Width Inside Handgrips.....	17" (431 mm)	17" (431 mm)
Overall Width.....	26" (665 mm)	26" (665 mm)
Overall Width Folded Upright.....	13" (330 mm)	13" (330 mm)
Overall Depth.....	28.5" (723 mm)	28.5" (723 mm)
Depth Folded Down Flat.....	36.4" (925 mm)	36.4" (925 mm)
Minimum Handle Height.....	33" (836 mm)	37" (935 mm)
Maximum Handle Height.....	38" (956 mm)	42" (1060 mm)
Seat to Floor Height.....	21.8" (556 mm)	26" (660 mm)
Approximate Height of Intended User.....	5'4" – 6'2"	6'2" – 6'8"
Maximum Weight Capacity		
EVENLY DISTRIBUTED.....	300 lb (136 kg)	300 lb (136 kg)
Product Weight.....	18.5 lb (8.4 kg)	19.1 lb (8.7 kg)
Seat Width.....	18" (460 mm)	18" (460 mm)
Seat Depth.....	9.4" (460 mm)	9.4" (460 mm)
Wheel Size.....	8" x 1.5" (201 mm x 38 mm)	8" x 1.5" (201 mm x 38 mm)