

Everest & Jennings[®]
by graham-field

WHEELCHAIRS



graham
field[™]



Everest & Jennings®

by graham-field

Value:

Graham-Field is known for a broad product offering at affordable pricing, without sacrifice to quality. Like you, we understand that the cost of a product is more than a price tag. The cost also includes the time and effort required to sell to and service your customers. At Graham-Field we are proud to offer the value found in the Everest & Jennings® line of wheelchairs.

Quality:

The E&J™ name has been synonymous with quality manual folding wheelchairs since 1933. In a time when many manufacturers are producing products at the lowest possible cost, often cutting corners on quality, Graham-Field is standing by its dedication to quality through product design, specification control, testing and inspection. It is easy to see why the E&J Traveler® name is synonymous with quality.

Service:

Our products are often only a day away with centralized distribution located in the metro areas of Atlanta, Georgia; St. Louis, Missouri; and Los Angeles, California. Graham-Field also carries a full line of E&J parts and upholstery to meet your service needs. In addition, our customer service department and technical service teams are always ready to help with your product needs.



Manual Folding Wheelchairs

Advantage® and Advantage® LX

3H01 & 3H02 Series

The Advantage, designed as a low-maintenance wheelchair, provides comfort and durability at an affordable price.

- Available in 16", 18" and 20" seat widths, 16" seat depth (additional 20"x18" on Traveler SE series)
- Black padded nylon and vinyl upholstery with standard chart pocket on back
- Carbon steel frame with hammertone finish provides strength and durability
- Tool-free adjustable front rigging with standard heel loops on swingaway footrests
- Comfortable desk and full length arms with padded armrests
- Maximum weight capacity 250 lb on Advantage series and 300 lb on Advantage LX series EVENLY DISTRIBUTED
- One (1) year limited warranty



K0001*

Traveler® SE

3E01 & 3E03 Series

The Traveler SE is the wheelchair of choice for standard transportation and adjustable seat height.

- Black vinyl upholstery with standard chart pocket on back
- Two-position axle and caster fork provide 2" seat height adjustment
- Maximum weight capacity 250 lb & 300 lb on 20"x18" models EVENLY DISTRIBUTED
- One (1) year limited warranty



K0001*
K0002
K0006

Traveler® LX

3E02 Series

The Traveler LX incorporates the legendary Everest & Jennings wheelchair chrome finish frame with saddle-joint tubular construction with dual axle. The premium blue vinyl upholstery is reinforced with an extra foam and heavy duty inner liner for enhanced comfort and durability.

- Chrome finish provides a long lasting and easy to clean finish
- Gray non-marking tires and aluminum footplates
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED
- Two (2) year limited warranty



K0001*
K0002

Bariatric Folding Wheelchairs

Traveler® HD

3G01 Series

The Traveler HD is a heavy duty wheelchair featuring a reinforced frame, double crossbrace, and a long lasting, embossed back and seat vinyl upholstery with a double inner liner for durability and support.

- Available in 20", 22" and 24" seat widths, 18" seat depth (1" longer than standard)
- Embossed black vinyl upholstery with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Reinforced frame at all critical load-bearing points
- Two-position axle and caster fork provides 2" seat height adjustment
- Tool-free adjustable front rigging with standard heel loops on footrests and durable aluminum black footplates
- Comfortable full and desk length arms with padded armrests
- Maximum weight capacity 500 lb EVENLY DISTRIBUTED
- Five (5) year limited warranty




K0007*

Traveler® HTC

3J01 Series

The Traveler HTC (Hospital Transport Chair) is an all-in-one solution wheelchair for the clinical transport need, as it incorporates a special attendant push handle and rear mounting attendant wheel locks for better caregiver assistance.

- Fixed desk arms and front rigging lock nuts ensure inventory control
- Anti-tippers and one black painted O2/IV telescopic holder for hospital/clinical settings come standard (black painted O2/IV pole is included after serial number HYT039WC)
- Maximum weight capacity 500 lb EVENLY DISTRIBUTED
- Five (5) year limited warranty



K0007*

Paramount XD

5PX10 Series

The Paramount XD, built for the bariatric wheelchair user, is constructed of 1" tubing to meet high weight capacity demands.

- Available in 26" and 30" seat widths, 19 ½" seat depth
- Comfortable desk arms with padded armrests
- Front rigging with extra wide aluminum footplate
- 8"x 2" non-marking solid casters and 24" solid rear wheels with chromed handrims
- Black reinforced nylon upholstery with standard chart pocket on back
- Maximum weight capacity of 650 lb EVENLY DISTRIBUTED
- Five (5) year limited warranty



K0007*

Specialty Folding Wheelchairs

Advantage® Recliner

3K01 & 3K02 Series

The Advantage Recliner, with a higher back and removable, full-support headrest for additional neck support, adjusts easily to fit every need.

- Available in 16", 18" and 20" seat widths, 17" seat depth and in 22" seat width, 19" seat depth
- Durable hammertone frame features manual full reclining back range from 90-170 degrees
- Comfortable desk and full length arms with padded armrests
- Embossed black vinyl upholstery with standard chart pocket on back
- Anti-tippers standard
- Metal spreader bar on back which works as anti-fold device
- Tool-free adjustable elevating front legrests
- Maximum weight capacity 300 lb EVENLY DISTRIBUTED on 16", 18" and 20" wide models (450 lb on 22" wide models)
- Five (5) year limited warranty



K0001* K0007*
E1226 E1226
E0966 E0966
E0971 E0971
K0195 K0195
E2201

Rehab Shower Commode

12022 & 12023 Series

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. Designed with a corrosion-resistant stainless steel frame, the chair can be positioned over a standard toilet or used as a portable, self-contained commode.

- Available with durable 24" rear wheels with handrims and solid rubber tires, with or without locking pin casters, or with four 5" casters. Wider models come with nylon tires for added durability
- Removable full length arms with composite armrests
- Available in Low-Back 14 ¼", 16" and Hi-Back 17 ¼", 19 ½"
- Available in 17" seat with 18" seat depth on Low-back 14 ¼" and 17 ¼" and the wider version 20" seat width on 16" Low-back and Hi-back 19 ½" models
- Sealed vinyl cushioned oval seat and frame side contours allow easy access and comfort in toileting activities and bowel management procedures
- Removable, hygienic plastic pail and lid are included
- Non-corrosive composite footplates
- Maximum weight capacity 250 lb EVENLY DISTRIBUTED on 14 ¼" and 17 ¼" back models (300 lb on 16" and 19 ½" back models)
- Two (2) year limited warranty



E0165*

Lightweight Folding Wheelchairs

Traveler® L3

3F01 Series

The Traveler L3 lightweight frame weighs less than 36 lb (without front rigging) making it the ideal wheelchair built for short and long term use.

- Available in 16", 18" and 20" seat widths, 16" seat depth
- Black nylon upholstery with standard chart pocket on back
- Non-marking polyurethane rear tires
- Hammertone finish and carbon steel frame provide strength and durability
- Comfortable padded armrests full and desk length
- Tool-free adjustable front rigging with standard heel loops on swingaway footrests
- Two-position axle and caster fork provide 2" seat height adjustment
- Maximum Weight Capacity 250 lb EVENLY DISTRIBUTED
- Five (5) year limited warranty



K0003*

Traveler® L3 Plus

3F012 Series

The Traveler L3 Plus is an improved version of the previous Traveler L3 series. The wheelchair comes standard with adjustable anti-tippers and a fully detachable 2" seat extension with hardware.

- Available in 16", 18", 20" seat widths including 2" removable seat extension kit converting the chair from 16" depth to 18" depth if necessary
- Includes fastener tab under seat connecting seat and back upholstery to eliminate seat-to-back gap
- Comfortable fixed flip-back desk length arms with padded armrests
- Maximum Weight Capacity 300 lb EVENLY DISTRIBUTED
- Five (5) year limited warranty



K0003*

Traveler® L4

3F02 & 3F03 Series

The Traveler L4 is a high-strength lightweight chair that weighs less than 34 lb (without front rigging) with an adjustable back.

- Available in 16", 18" and 20" seat widths, 16" seat depth. Also 18" and 20" seat widths, 18" seat depth
- Quick release axle and height adjustable desk arms are available for some models
- Black nylon upholstery with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Anti-tippers come standard
- Comfortable flip-back desk and full length arms with padded armrests
- Maximum weight capacity 250 lb EVENLY DISTRIBUTED (300 lb on models with 18" seat depth)
- Five (5) year limited warranty



K0004*

Transport Chairs

Aluminum Transport Chair

EJ78X-1 & EJ79X-1 Series

The Everest & Jennings Transport Chair is designed as an affordable short-term patient transfer. Its compact size and lighter weight allow a companion to lift and store the chair with ease.

- 17" and 19" seat widths (Deluxe model only available in 19")
- Lightweight aluminum frame available in Aqua, Black, Red, Blue, Silver and Pink (Black and Blue Celeste are available for Deluxe Models only)
- Fold-down back makes transport and storage easy
- Black, padded, nylon upholstery is flame retardant
- Tool-free adjustable detachable swingaway footrests
- Comfortable fixed full length arms with padded armrests
- Comes standard with 3" leg adjustable loops and strap to slide over installed footrests
- Includes carry straps on seat for lifting and folding ease
- 42" length positioning belt
- Rear mounted attendant locks on rear casters make it easy for a companion to engage
- Weighs 20 lb (without front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- One (1) year limited warranty

colorSelect®

E1038*



Steel Transport Chair

EJ79X-1 Series

- Steel frame is finished in a durable hammertone coating
- Weighs 24 lb (without front rigging)
- Maximum Weight Capacity 250 lb EVENLY DISTRIBUTED
- One (1) year limited warranty



E1038*

Deluxe Transport Chair

EJ87X-1 Series

- Handle mounted attendant locks
- 12" flat-free rear wheels with 8" swivel front casters
- Weighs 21.38 lb (without front rigging)
- Maximum Weight Capacity 300 lb EVENLY DISTRIBUTED
- One (1) year limited warranty



Bariatric Transport

This chair satisfies additional compact transport needs with extra capacity. Handle-mounted attendant wheel locks are easily engaged by a companion.

- Dual crossbrace on gusseted aluminum frame with silver epoxy painted finish
- Available in 24" width, 18" seat depth
- 12" flat-free rear wheels with 8" swivel front casters
- Removable flip-back desk length padded arms provide user comfort
- Detachable swingaway footrests
- Available with anti-tippers on models after serial number IOS291
- Weighs 31 lb (without front rigging)
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty
- Item No. EJ777-3




E1039*

Accessories

Anti-tippers	
	
90763565	Anti-tipper, Advantage
90763065	Anti-tipper, Traveler SE & HD
90763065	Anti-tipper, Traveler L3
90763230A	Anti-tipper, Traveler L3 Plus

IV/O2 Holders	
	
B10320F	O2/IV Holder Combination (Advantage, SE, L3, L3 Plus, L4, HD, LX, Paramount XD)
B10325	O2 Holder (Advantage, SE, L3, L3 Plus, L4, HD, LX, Paramount XD)


Positioning Belts	
	
92133040	Positioning Belt, Metal Buckle
300-505	Positioning Belt, Auto Style
EJ76-BELT	Positioning Belt, Plastic Buckle

O2 Cylinder Holder Bag	
	
EJ02BAG07	O2 Bag

Footrest	
	
90763030	Footrest with Composite Footplate
90763430	Footrest with Aluminum Footplate

Legrest	
	
90763031	Elevating Legrest with Composite Footplate
90763431	Elevating Legrest with Aluminum Footplate

Articulating Legrest	
	
90763432	Heavy duty articulating legrest complete w/ aluminum footplate & panel

Padded Seat	
	
90835070	Oval cutout seat for 17 1/4" and 14 1/4"
90835071	Oval cutout seat for 16" and 19"
90835175	Open front (C-shaped) for 17 1/4" and 14 1/4"

Item No.	Description
Advantage	
3H010120 / 3H010130	Advantage 18X16 Desk Arms, Footrest / Legrest
3H020120 / 3H020130	Advantage LX 20X16 Desk Arms, Footrest / Legrest
Traveler SE	
3E010220 / 3E010230	Traveler SE 16X16 Desk Arms, Footrest / Legrest
3E010120 / 3E010130	Traveler SE 18X16 Desk Arms, Footrest / Legrest
3E010320 / 3E010330	Traveler SE 20X16 Desk Arms, Footrest / Legrest
Traveler LX	
3E020120 / 3E020130	Traveler LX 18X16 Desk Arms, Footrest / Legrest
Traveler L3	
3F010120 / 3F010130	Traveler L3 18X16 Desk Arms, Footrest / Legrest
3F010320 / 3F010330	Traveler L3 20X16 Desk Arms, Footrest / Legrest
Traveler L3 Plus	
3F012220 / 3F012230	Traveler L3 PLUS 16X16 Desk Arms, Footrest / Legrest
3F012120 / 3F012130	Traveler L3 PLUS 18X16 Desk Arms, Footrest / Legrest
3F012320 / 3F012330	Traveler L3 PLUS 20X16 Desk Arms, Footrest / Legrest
Traveler L4 (K0004)	
3F020220 / 3F020230	Traveler L4 16X16 Desk Arms, Footrest / Legrest
3F020120 / 3F020130	Traveler L4 18X16 Desk Arms, Footrest / Legrest
3F020320 / 3F020330	Traveler L4 20X16 Desk Arms, Footrest / Legrest
Traveler HD Heavy Duty	
3G010420 / 3G010430	Traveler HD 22X18 Desk Arms, Footrest / Legrest
3G010520 / 3G010530	Traveler HD 24X18 Desk Arms, Footrest / Legrest
Traveler Hospital Transport Chair (HTC)	
3J010330	Traveler HTC 20X18 Desk Arms, Elevating Legrest
3J010430	Traveler HTC 22X18 Desk Arms, Elevating Legrest
3J010530	Traveler HTC 24X18 Desk Arms, Elevating Legrest
Paramount XD Heavy Duty	
5PX10620 / 5PX10630	Paramount XD 26W Desk Arms, Footrest / Legrest
5PX10820 / 5PX10830	Paramount XD 30W Desk Arms, Footrest / Legrest
Aluminum Transport Chairs	
EJ781-1 / EJ782-1	Transport Chair Red 17" / 19"
EJ785-1 / EJ786-1	Transport Chair Black 17" / 19"
EJ787-1 / EJ788-1	Transport Chair Blue 17" / 19"
EJ790-1 / EJ791-1	Transport Chair Pink 17" / 19"
EJ793-1	Transport Chair Aqua 19"
EJ870-1 / EJ871-1	Deluxe Transport Chair Black / Blue Celeste
Steel Transport Chair	
EJ796-1	Transport Chair Silver Vein Steel 19"
Bariatric Transport	
EJ777-3	Transport Bariatric Chair 24W
Advantage Reclining Wheelchair	
3K010230 / 3K010250	Advantage Recliner 16X17 Legrest, Desk Arms / Full Arms
3K010130 / 3K010150	Advantage Recliner 18X17 Legrest, Desk Arms / Full Arms
3K010330 / 3K010350	Advantage Recliner 20X17 Legrest, Desk Arms / Full Arms
3K020430 / 3K020450	Advantage Recliner 22x19 Legrest, Desk Arms / Full Arms
Rehab Shower Commode	
12022010 / 12022011	14 1/4" LoBack / 17 1/4" LoBack, Caster Pin Locks
12022013	14 1/4" LoBack with 5" Rear Locking Caster
12022110 / 12022111	17 1/4" HiBack / 17 1/4" HiBack, Caster Pin Locks
12023010 / 12023110	16" LoBack / 19" HiBack

*For a complete list of sizes and options available, please visit www.grahamfield.com

Follow us on   

1.770.368.4700

Information contained herein is subject to change. The most current and complete product information can be found on our website.

www.grahamfield.com

© 2014, GF Health Products, Inc. All Rights Reserved
 Advantage, E&J, Everest & Jennings, Graham-Field, and Traveler are trademarks of GF Health Products, Inc. GF Health Products, Inc. is an ISO 13485 Certified Company.