

**The John Bunn Neb-u-Lite EV2 Nebulizer comes complete with:**

- 1 Compressor
- 1 disposable nebulizer cup (reusable for single patient use only)
- 1 Mouthpiece
- 1 PVC tubing
- 5 Spare filters
- 1 Adult and 1 Pediatric Mask

\* Boiling nebulizer components is not recommended



**Features:**

- Small footprint and light weight nebulizer compressor at a competitive price point for providers
- Compressor weighs only 2.7lb (1.2 kg)
- Reliable, piston-driven compressor is efficient and delivers consistent particle size of <5 microns
- Built-in nebulizer cup holder for hands-free use
- Complete with disposable nebulizer kit
- 5 year limited warranty
- Powerful piston pump is efficient, dependable and quiet
- Designed to meet the requirements of HCPCS Code: E0570



Item	UOM
JB0112-090	10 ea/cs

\* Federal Law (USA) restricts this device to sale by or on the order of a physician. This product should not be used unless the operator has been instructed by a qualified health care professional.

**Technical Specifications**

Method of Nebulization .....	Piston Driven
Power .....	120VAC, 60Hz
Operating Temperature Range .....	50°F to 104°F (10°C to 40°C)
Operating Humidity Range .....	30% to 95% RH
Storage Temperature Range .....	-40°F to 158°F (-40°C to 70°C)
Storage Humidity Range .....	10 to 95% RH
Weight .....	2.7 lb (1.2 kg)
Width .....	5.8"
Height .....	3.3"
Depth .....	5.5"
Max Pressure (Compressor) .....	30 PSI
Flow Range (Compressor) .....	0 - 8 liters/min
Particle Size (Nebulizer) .....	0.5 - 5 microns
Capacity (Nebulizer) .....	8 ml