

Labtron's Rechargeable Otoscope Ophthalmic Set offers a Double-Seat Charger Base designed to hold and charge the included Otoscope and Ophthalmoscope. Each scope features a durable, ergonomic handle equipped with an LCD screen indicating brightness level and battery life.

### Standard Features:

- Double-Seat Charger Base with LED indicator to show charge
- 2 Durable, Ergonomic Handles with LCD screen
- LED bulb with dimmable rheostat for brightness control
- Otoscope Head with extra tips
- Ophthalmoscope Head
- 5VDC Power Adapted and USB Cable
- Insufflation Bulb
- 4 rechargeable 1600mAh Batteries
- One (1) Year Limited Warranty



Item	Description	Shipping Wt.	UOM
1227-3	Rechargeable Otoscope Ophthalmic Set	2.5 lb	1 ea

