

Designed to provide assistance to those having difficulty bathing.



Features:

- Comfortable, cushioned seat and backrest are water tight and easy to clean
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Knock-down design makes transporting easy
- Quick and easy tool-free assembly
- Maximum Weight Capacity
300 lb - EVENLY DISTRIBUTED



Item#	Description	UOM
7944KD-1	Padded Bath Seat	1 EA
7945KD-1	Padded Bath Seat with Arms	1 EA
7946KD-1	Padded Commode Bath Seat	1 EA
7947KD-1	Padded Commode Bath Seat with Arms	1 EA

Specifications:

Overall Width.....	20.5"
Overall Depth.....	21"
Seat Width.....	15.75"
Seat Depth.....	16"
Minimum Seat Height.....	18"
Maximum Seat Height.....	22"
Maximum Weight Capacity	300 lb EVENLY DISTRIBUTED