

Convenient user-friendly snap-on design provides easy access for drink bottle



Standard Features:

- Easily snaps and adjusts by tension clip-ons onto walkers, rollators, wheelchairs, transport chairs, quad and standard canes with tubing diameter from 5/8" to 1" (1.58 cm to 2.54 cm)
- Tension adjustable clamp bottle holder united to a 3 1/2" (8.89 cm) diameter shelf base by two screws securing two metal clip-ons for extra support
- Frame made with durable high-density polypropylene
- 3 1/2" (8.89 cm) diameter composite shelf base to support bottle
- Total product weight: 2.6 ounces (73.70 g)
- Product dimensions: 3" x 3 1/2" x 3" (7.62 cm x 8.89 cm x 7.62 cm)
- Fits 2 1/2"-3" (6.35 cm x 7.62 cm) in diameter bottles (holder can be adjusted to fit bottle size)
- Weight capacity: bottled beverage containers from 8 to 20 oz (236.59 ml to 591.47 ml) evenly distributed – Maximum diameter 2.5" to 3" (6.35 cm x 7.62 cm)



Warranty:

Six (6) months limited warranty

Item	Description	UOM
EJCUPH19	Cup Holder	1 ea



**This product is not made with natural rubber latex*