

EXAM ROOM SUPPLIES

STEP STOOLS

Features

- Non-skid, ribbed rubber mat platform
- Rubber tipped legs
- One (1) year limited warranty



GF1841-2



GF1840-2



GF1842-1

Specifications:

		overall Height		Height		Width		Depth		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
GF1840C-2	Step-Up Stool Chrome Plated Steel	35.75	90.8	9"	22.9	14"	35.6	11"	27.9	250	113.4	2/CS
GF1841-2	Heavy-Duty Foot Stool Chrome Plated Steel			9"	22.9	14"	35.6	11"	27.9	350	158.8	2/CS
GF1840-2	Economy Foot Stool Epoxy Finished Steel			9"	22.9	14"	35.6	11"	27.9	250	113.4	2/CS
GF1842-1	Bariatric Foot Stool			9"	22.9	14"	35.6	11"	27.9	500	226.8	1 ea

WARNINGS: Cancer and Reproductive Harm -
www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.
www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Grafco is a trademark of GF Health Products, Inc.
GF Health Products, Inc. is an ISO 13485 Certified Company.