



RAISED TOILET SEATS & SAFETY RAILS

RAISED TOILET SEAT CUSHION

Essential item to assist with sitting and standing from a toilet seat.

Features

- Adds 4" to toilet seat height making sitting and standing easier
- Soft and comfortable foam with durable cover
- 4 straps securely hold cushion to the toilet seat



Specifications:

	overall Length		overall Depth		overall Height		UOM
	inch	cm	inch	cm	inch	cm	
6484R-2	13.4	34	15.75	40	4	10.1	1 ea

⚠️ WARNINGS: Cancer and Reproductive Harm -
www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Lumex is a trademark of GF Health Products, Inc.
GF Health Products, Inc. is an ISO 13485 Certified Company.