

## OXIMETERS

# FINGER PULSE OXIMETER

The JB02020 Finger Pulse Oximeter provides simple, fast, and accurate oxygen saturation (SpO<sub>2</sub>) and pulse rate measurements. Its lightweight and compact design makes it an ideal solution for home, clinical and EMS settings.

### Features

- Lightweight, pocket-size solution
- Provides fast and accurate SpO<sub>2</sub>, pulse rate and Perfusion Index measurements
- Large LED digital display that is easy to read
- Simple one button operation
- Automatic power off function to conserve batteries
- Accommodates widest range of finger sizes from pediatric to adult

### Specifications:

**Display:** Oxygen saturation (SpO<sub>2</sub>), Pulse rate and Perfusion Index

**Display Range:** SpO<sub>2</sub>: 35% to 100%

**Pulse Rate:** 30 to 240 bpm

**Accuracy (% SpO<sub>2</sub>):** 80% - 99%: ±2% 70% - 100%: ±2%

**Pulse Rate Accuracy:** ±2%

**Perfusion Index:** 0-20%

**Alarms:**

SpO<sub>2</sub> low limit: 90%

Pulse Rate:

High limit: 120 bpm

Low limit: 50 bpm

Low Battery indicator

**Unit Weight:** .1 lb (52g)



	overall Length		overall Width		overall Height		UOM
	inch	mm	inch	mm	inch	mm	
JB02020	2.2	56	1.3	34	1.2	30	1 ea

**WARNINGS:** Cancer and Reproductive Harm - [www.P6gWarnings.ca.gov](http://www.P6gWarnings.ca.gov)

Color variations may occur due to monitor or printer settings.

GF Health Products, Inc. is an  
ISO 13485:2016  
MDSAP Certified Company  
**MDSAP**  
MEDICAL DEVICE SINGLE AUDIT PROGRAM

© 2020, GF Health Products, Inc. All Rights Reserved.  
Hausted is a registered trademark of GF Health Products, Inc.



© 2020, GF Health Products, Inc. All Rights Reserved.  
John Bunn is a trademark of GF Health Products, Inc.  
GF Health Products, Inc. is an ISO 13485 Certified Company.

**graham field**

GF200016RevA20