



HEMOCARE BED ACCESSORIES

HALF LENGTH CLAMP-ON BED RAIL

The Half Clamp-On Bed Rail provides a full head coverage of 32 inches, designed to offer bed protection and eliminating gaps in between the railing. The non-telescoping tube structures offers a very sturdy rail that can be used as a secondary handle and support.



Features

- Improved rail design decreases the spacing between the railing thus reducing the risk of entrapment
- Tool-free installation so it easily mounts to most bed frames using convenient clamp (not for use with "Pan-Style" bed surfaces)
- Intended to define the edge of the sleep surface and to help prevent an individual from inadvertently rolling out of bed
- Designed to be used with the Lumex Patriot™ beds that are commonly used as standard-height or low beds
- Heavy-duty steel construction
- Rails are powder-coated with an environmentally-friendly, long-lasting finish and secure gripping surface
- 3 Year Limited Warranty
- Designed to meet the requirements of HCPCS E0305



Specifications:

	overall Length		overall Height		product Weight		UOM
	inch	mm	inch	mm	lb	kg	
GF6650B-1	32	813	18	457	16.8	7.6	1 pr

WARNINGS: Cancer and Reproductive Harm - www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



© 2020, GF Health Products, Inc. All Rights Reserved.
Lumex is a trademark of GF Health Products, Inc.
GF Health Products, Inc. is an ISO 13485 Certified Company.