



RAISED TOILET SEATS

DELUXE TOILET SEAT RISERS

Essential item to assist with sitting and standing from a toilet seat.

**Features**

- Adds 3.5" to toilet seat height, reducing bending effort
- Securely attaches to toilet bowl using bolts to prevent shifting
- Designed to work with your existing toilet seat and lid to better blend in with decor
- Available in round or elongated toilet bowl sizes

**With Arms (6482R-2, 6482ER-2)**

- Push buttons allow for armrest to be added or removed easily

**Flip-up Seat (6481R-2, 6481ER-2)**

- Designed to allow toilet seat riser to flip-up for easy cleaning



**Specifications:**

	overall Width		Seat Depth		Seat Height		Width Between Arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6482R-2	22.5	57.15	17.25	43.82	3.5	8.89	19.5	49.53	300	136.08	2/cs
6482ER-2	23.5	59.69	19	48.26	3.5	8.89	20.5	52.07	300	136.08	2/cs
6481R-2	13.8	35.05	17.5	44.45	3.5	8.89	N/A	N/A	300	136.08	2/cs
6481ER-2	13.8	35.05	19	48.26	3.5	8.89	N/A	N/A	300	136.08	2/cs

**WARNINGS:** Cancer and Reproductive Harm - [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Color variations may occur due to monitor or printer settings.



© 2020, GF Health Products, Inc. All Rights Reserved.  
Hausted is a registered trademark of GF Health Products, Inc.



© 2020, GF Health Products, Inc. All Rights Reserved.  
Lumex is a trademark of GF Health Products, Inc.  
GF Health Products, Inc. is an ISO 13485 Certified Company.