



OVERBED TABLES

Manufactured with the lowest base in the industry, our Overbed tables move easily under most hospital beds or recliners.

Features

- Lowest profile base on the market
- Exposed steel base with no plastic cover allows for easy cleaning of base and casters
- Laminate or Thermofoil top that provides durability
- Available in 4 attractive finishes: Candlelight, Hardrock Maple, Medium Cherry, and Medium Oak
- Easy, one-hand activated secondary surface reduces infection control issues and eliminates hidden spaces between surfaces for bacteria to collect
- Integrated cup holders
- Optional Vanity drawer with mirror slides out from either side of the table
- Integrated accessory rail for urinals, storage baskets etc.
- Shock-absorbent gas spring cylinder for easy height adjustment
- Two (2) Year Limited Warranty



Specifications:

| | | table top Length | | table top Width | | Height Range | | Product Weight | | Weight Capacity EVENLY DISTRIBUTED | | UOM |
|--------------|--------------------------------------------------------------|-------------------|-------------------|-----------------|-------|--------------|----------|----------------|-------|---------------------------------------|-------|------|
| | | inch | cm | inch | cm | inch | cm | lb | Kg | lb | Kg | |
| GF9000-LSTMO | Single Top, Laminate, Medium Oak, Graphite Grey U Base | 32.5 | 82.55 | 19 | 48.26 | 28-44 | 71-111.8 | 34 | 15.42 | 100 | 45.36 | 1/ea |
| GF9000-LSTMC | Single Top, Laminate, Medium Cherry Graphite Grey U Base | 32.5 | 82.55 | 19 | 48.26 | 28-44 | 71-111.8 | 34 | 15.42 | 100 | 45.36 | 1/ea |
| GF9000-LSTHM | Single Top, Laminate, Hardrock Maple, Graphite Grey U Base | 32.5 | 82.55 | 19 | 48.26 | 28-44 | 71-111.8 | 34 | 15.42 | 100 | 45.36 | 1/ea |
| GF9000-LSTCL | Single Top, Laminate, Candlelight, Graphite Grey U Base | 32.5 | 82.55 | 19 | 48.26 | 28-44 | 71-111.8 | 34 | 15.42 | 100 | 45.36 | 1/ea |
| GF9000-TSTMO | Single Top, Thermofoil, Medium Oak, Graphite Grey U-Base | 32.5 | 82.55 | 19 | 48.26 | 28-44 | 71-111.8 | 34 | 15.42 | 100 | 45.36 | 1/ea |
| GF9000-TSTMC | Single Top, Thermofoil, Medium Cherry Graphite Grey U-Base | 32.5 | 82.55 | 19 | 48.26 | 28-44 | 71-111.8 | 34 | 15.42 | 100 | 45.36 | 1/ea |
| GF9000-TSTHM | Single Top, Thermofoil, Hardrock Maple, Graphite Grey U-Base | 32.5 | 82.55 | 19 | 48.26 | 28-44 | 71-111.8 | 34 | 15.42 | 100 | 45.36 | 1/ea |
| GF9000-TSTCL | Single Top, Thermofoil, Candlelight, Graphite Grey U-Base | 32.5 | 82.55 | 19 | 48.26 | 28-44 | 71-111.8 | 34 | 15.42 | 100 | 45.36 | 1/ea |
| GF9000-TDTMO | Flip Top, Thermofoil, Medium Oak, Graphite Grey U-Base | 32+13.5 (flip)=46 | 82.55+34.3=116.84 | 19 | 48.26 | 28-44 | 71-111.8 | 41 | 18.6 | 100 | 45.36 | 1/ea |
| GF9000-TDTMC | Flip Top, Thermofoil, Medium Cherry Graphite Grey U-Base | 32+13.5 (flip)=46 | 82.55+34.3=116.84 | 19 | 48.26 | 28-44 | 71-111.8 | 41 | 18.6 | 100 | 45.36 | 1/ea |
| GF9000-TDTHM | Flip Top, Thermofoil, Hardrock Maple, Graphite Grey U-Base | 32+13.5 (flip)=46 | 82.55+34.3=116.84 | 19 | 48.26 | 28-44 | 71-111.8 | 41 | 18.6 | 100 | 45.36 | 1/ea |
| GF9000-TDTCL | Flip Top, Thermofoil, Candlelight, Graphite Grey U-Base | 32+13.5 (flip)=46 | 82.55+34.3=116.84 | 19 | 48.26 | 28-44 | 71-111.8 | 41 | 18.6 | 100 | 45.36 | 1/ea |
| GF9000-BSK | Storage Basket - 10X6x9 With Adapter | | | | | | | | | | | 1/ea |
| GF9000-VNTY | Vanity Drawer With Mirror | | | | | | | | | | | 1/ea |
| GF9000-UH | Urinal Hook With Adapter | | | | | | | | | | | 1/ea |
| GF9000-CST | Non-Locking Caster | | | | | | | | | | | 1/ea |

WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com