

SPHYGMOMANOMETERS

# AUTOMATIC BLOOD PRESSURE MONITOR

Blood pressure management is one of the most important things you can do to maintain your health and is essential to preventing heart disease and strokes. More doctors are recommending the use of home blood pressure monitors so that patients can regularly, conveniently check their blood pressure.

Lumiscop's Automatic Upper Arm Blood Pressure Monitor makes home blood pressure management stress free with its ease of use and portable design.

**Features**

- Automatic Inflation and Deflation Blood Pressure Monitor
- Digital LCD
- Memory feature - 2 user with 90 sets of memories each
- Average value of the latest 3 readings
- Uses 4 AA batteries (not included)
- Includes D-ring cuff designed to be self-applied
- Irregular heartbeat detection



**Specifications:**

	overall Width		overall Length		overall Height		UOM	
	inch	mm	inch	mm	inch	mm		
	11.33	4.72	120	4.13	105	2.17	55	1/ea

**WARNINGS:** Cancer and Reproductive Harm - [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Color variations may occur due to monitor or printer settings.