

PHLEBOTOMY CHAIRS



PHLEBOTOMY CHAIRS

Hausted phlebotomy chairs were designed with both the patient and phlebotomist in mind. The sensible and affordable design offers a small footprint and modern style while keeping patient comfort and well-being the highest priority.

- Powder coat frame finish - Foghorn
- Fully upholstered padded arm rests
- One flip arm standard (Dual flip arm option on 6200 series)
- Height adjustable arm rests



6200 Series Bariatric

6200-xx upholstered bariatric blood draw/phlebotomy chair with padded arm rest. Dual Flip Arm shown on cover



6100

6100-xx upholstered blood draw/phlebotomy chair with padded arm rest

Standard Color Options (-xx)

Available in 24 color choices to complement most decors.

| | | | | | |
|-------------------|------------------|--------------|-----------------|-------------------|------------------|
| WH Super White | DO Doe | PM Palm | SH Shadow | AR Artichoke | TE Teal |
| NO Northwoods | DB Dusty Blue | PU Plum | WW Wedgewood | BT British Tan | GR Gray |
| AL Alabaster | ST Saddle Tan | PL Pearl | TL Tea Leaf | MO Moss | OG Ocean Gray |
| AQ Aqua Green | BU Blue | VI Violet | NA Navy | WG Warm Gray | BL Black |

Swatch Card is available from your representative

Specifications:

| | 6100 | | 6200 | |
|--|-----------|--------------|-----------|--------------|
| | inch | cm | inch | cm |
| Overall Width: | 28.50 | 72.39 | 34.00 | 86.36 |
| Overall Depth: | 29.00 | 73.66 | 29.00 | 73.66 |
| Overall Height: | 37.50 | 95.25 | 37.50 | 95.25 |
| Seat Width: | 18.50 | 46.99 | 24.50 | 62.23 |
| Seat Depth: | 18.50 | 46.99 | 18.50 | 46.99 |
| Seat Height: | 22.00 | 55.88 | 22.00 | 64.77 |
| Armrest Height: -low adjustment | 27.00 | 62.23 | 27.00 | 62.23 |
| Armrest Height: -high adjustment | 36.00 | 66.04 | 36.00 | 66.04 |
| Back Width: | 19.50 | 49.53 | 25.50 | 64.77 |
| Back Height: | 17.00 | 43.18 | 17.00 | 43.18 |
| Weight Capacity: EVENLY DISTRIBUTED | lb 600 | kg 272.16 | lb 800 | kg 362.87 |
| Shipping Weight: | 55 | 24.95 | 55 | 24.95 |

| Options: available at additional cost | item number | item number |
|--|-------------|-------------|
| | None: | 6100-xx |
| Dual Flip Arms: | — | 6200-DA-xx |

— Denotes Feature Not Available / Does Not Apply

WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP) which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to - www.P65Warnings.ca.gov/furniture

Color variations may occur due to monitor or printer settings.

GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company
MDSAP
Intertek



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com

Product Hyperlinks Enabled