



SUPPORT SURFACES


SELECT FOAM MATTRESSES

LS100 Series

- High Density Foam Mattress
- Convoluted Topper for Excellent Airflow
- Durable Nylon/Vinyl Cover
- Compress rolled for cost effective storage and freight
- ★ See Chart below for Links, Item Numbers, Specifications, Options and Warranty

LS150 Series

- High Density Foam Mattress
- Firm Side Rails and Soft Heel Section
- Long Term Care Replacement Mattress
- Stretch Urethane/Vinyl Cover
- Compress rolled for cost effective storage and freight
- Available with Raised Perimeter Bolster (LS150-35P and LS150-42P)
- ★ See Chart below for Links, Item Numbers, Specifications, Options and Warranty



Raised Perimeter Bolster
for patients likely to fall out of or exit the bed unattended

Firm Side Rails
aids in fall prevention during ingress / egress

Convoluted Topper
creates excellent airflow

Soft Heel Section
for patients with sensitive feet

Specifications

	Width		Length		Height		Weight Capacity EVENLY DISTRIBUTED		Durable Nylon/Vinyl Stretch Urethane/Vinyl		High Density Foam Convoluted Topper		Firm Side Rails	Soft Heel Section	Raised Perimeter	Pressure Redistribution	Meets 16 CFR 1633	Mattress Warranty in Months	Cover Warranty in Months
	inch	cm	inch	cm	inch	cm	lb	kg	cover	build									
LS100-35	35.00	88.90	80.00	203.20	6.00	15.24	350	158.76	■	—	■	■	—	—	—	■	■	60	12
LS100-42	42.00	106.68	80.00	203.20	6.00	15.24	350	158.76	■	—	■	■	—	—	—	■	■	60	12
LS150-35	35.00	88.90	80.00	203.20	6.00	15.24	400	181.44	—	■	■	—	■	■	—	■	■	60	12
LS150-42	42.00	106.68	80.00	203.20	6.00	15.24	400	181.44	—	■	■	—	■	■	—	■	■	60	12
LS150-35P	35.00	88.90	80.00	203.20	6.00	15.24	400	181.44	—	■	■	—	■	■	■	■	■	60	12
LS150-42P	42.00	106.68	80.00	203.20	6.00	15.24	400	181.44	—	■	■	—	■	■	■	■	■	60	12

■ Denotes Standard Feature

— Denotes Feature Not Available/Does Not Apply



Shipping Box Dimensions:
40 x 14 x 14 inches
(101.60 x 35.56 x 35.56 cm)

Color variations may occur due to monitor or printer settings.

WARNING: Cancer and Reproductive Harm - www.P66Warnings.ca.gov



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

