



PAIN MANAGEMENT

TENS is a device that uses electrical current for stimulating nerves to help relieve pain. Usually used as an alternative to medications. EMS is an electronic device used to deliver electrical impulses to the body in order to make the muscles contract.



GF-3

Analog TENS

Prescription Required



- Feature:** Dual Channel
- Modes:** 3 - Standard, Burst, Modulation
- Includes:** TENS Unit, soft carry case, instructions, four electrodes and lead wires
- Uses:** 9 Volt Battery (NOT INCLUDED)

GF-3T

Analog TENS with Timer

Prescription Required



- Features:** Dual Channel, Timer
- Modes:** 3 - Standard, Burst, Modulation
- Includes:** TENS Unit, soft carry case, instructions, four electrodes and lead wires
- Uses:** 9 Volt Battery (NOT INCLUDED)

GF-TX5

Digital TENS

Prescription Required



- Features:** Dual Channel, Timer, LCD Screen
- Modes:** 3 - Standard, Burst, Modulation
- Includes:** TENS Unit, soft carry case, instructions, four electrodes and lead wires
- Uses:** 9 Volt Battery (NOT INCLUDED)

GF-DF5

Digital TENS/EMS

Prescription Required



- Features:** Dual Channel, Dual Function TENS/EMS, Timer, LCD Screen
- Modes:** 3 - Standard, Burst, Modulation
- Includes:** TENS Unit, soft carry case, instructions, four electrodes and lead wires
- Uses:** 3 AAA Batteries (NOT INCLUDED)

WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov

GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company
MDSAP
Intertek MEDICAL DEVICE SINGLE AUDIT PROGRAM

© 2021, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

Product Hyperlinks Enabled

Color variations may occur due to monitor or printer settings